

## Calibration Procedure:

1. Start with the wristwatch secured on the wrist, insert a micro-SD card into the slot, and turn the switch on.
2. Keep the wrist in a horizontal position with the watch facing up.
3. Rotate wrist 180 degrees CCW
4. Rotate wrist 180 degrees CW
5. Repeat 2 + 3 x2 with 1 second movements
6. Repeat 2 + 3 x2 with 2 second movements

1. Start with elbow bent, the arm facing upwards, and the wrist facing left
2. Rotate arm around elbow 90 CW, towards the right arm
3. Rotate arm around elbow 90 CCW, away from the right arm
4. Repeat 2 + 3 x2 with 1 second movements
5. Repeat 2 + 3 x2 with 2 second movements

1. Start with the arm and wrist in a horizontal position with the watch facing up.
2. Rotate arm around elbow 90 CW, towards the chest, maintaining a horizontal position.
3. Rotate arm around elbow 90 CCW, away from the chest
4. Repeat 2 + 3 x2 with 1 second movements
5. Repeat 2 + 3 x2 with 2 second movements

1. Start with the wrist facing to the left and the arm stretched away from the body to the left.
2. Rotate arm around shoulder horizontally 90 CW, towards the right arm.
3. Rotate arm around shoulder horizontally 90 CCW, away from the right arm
4. Repeat 2 + 3 x2 with 1 second movements
5. Repeat 2 + 3 x2 with 2 second movements

Turn the switch off, take the wristwatch off, and take the micro-SD card out of the device.